



WWW.SENSICA.COM

MD 04/18/ A3

SENSICA



Sensilift

User Guide | Mode D'emploi | Benutzerhandbuch
Istruzioni D'utilizzo | Instrucciones De Uso

01 Getting Started

Hi, welcome to Sensica. We're so excited for you to start your journey to younger looking skin. Before you do your first treatment, check out this Quick Guide with extra tips for best results.

1. 48 hours before starting your first full treatment, perform a preliminary test on each area you intend to treat.
2. Plug the device into an electric socket. You should see the "status indicator" blinking yellow. (front light)
3. Make sure your skin is clean and dry, free from creams and make up.
4. Apply an even layer of Base Gel on the area intended for treatment.
5. Set the energy level by pressing the RF button (see below) repeatedly, until you reach the desired level: 1, 2 or 3.
6. Press the Massage button (see below) to set the desired massage level: 1 or 2.

RF Button - ☺

Massage Button - ☺

Hey ! I'm Sensica's beauty consultant, and I'm here to help you with some extra tips...

**Remember, patience is key
Always do a preliminary test before beginning a full treatment.
Start at RF level 3, and reduce the energy level to suit your comfort.**

The massage feature is all about getting into the deep layers of the skin - it's what makes this device so effective.

02 Treatment Schedule

1. Place the RF electrodes (massage fingers) in full contact with your skin.
2. You should see a steady green light in the "status indicator" as soon as there is full contact with the skin.
3. Move the device over each treatment area in circular motions for 5 minutes (10 minutes for forehead).
4. The treatment will create a warm sensation on your skin.
5. Wipe off any Base Gel residue from the treated area. To enhance results, apply Sensica "Age Defying" cream (also available at www.sensica.com).
6. For best results, use Sensilift once a week on each area for the period of 8 weeks.

You can use a mirror during treatment to keep an eye on the light.

Don't be tempted to treat your whole face in one go - stick to areas the size of a business card.

The safety mechanism in the device will regulate the RF energy level to suit your skin temperature, so don't worry if you feel the heat come and go throughout the treatment.

**That's it!
Your skin will feel lifted,
smoother and tighter.
Enjoy!**

03 Optimizing Results

1. Results of Sensilift treatments may differ between individuals, depending on environmental and biological factors, heredity and age.
2. The device delivers RF energy only when the green light is steady. If this doesn't happen after a few seconds, try these steps:
 - make sure you have selected an RF level.
 - make sure the RF electrodes (massage fingers) are in full contact with the skin.
 - make sure you have applied enough Base Gel to facilitate good contact.
3. If you don't feel heat being delivered to the treatment area during treatment:
 - make sure the status indicator has a steady green light.
 - increase the RF level.
 - decrease the treatment area and reduce the radius of your circular motion.
4. To sustain optimal long term results, maintenance treatments might be performed every 4-8 weeks, as needed.

Remodelling and changes in collagen production may take time. Be patient - you will see results within 8 treatments.

The secret to getting the best result is to stick with the protocol and persevere.

Feel free to ask me any question about your skin or this treatment at: sensicabeautyconsultant@sensica.com

Warranty Card

Congratulations on your Sensica purchase!

We have every faith in the quality of each Sensica device – we know they meet the highest standards, because we subject them to rigorous quality and performance testing. In the unlikely event that you do find a defect in a Sensica device, it may be covered by this warranty.

Two Year Limited Warranty

This warranty covers only Sensica devices. Other Sensica products, such as pads or creams, are defined as consumables and therefore excluded from warranty coverage.

This limited warranty covers manufacturer defect related problems in your Sensica device, if you purchased it directly from www.sensica.com or from an authorized partner. Warranty approvals will only be granted

where an original receipt is provided. This warranty is limited to devices proven to be defective.

In order to obtain warranty service on your Sensica device, please go to our website or call our customer service center for assistance (you will find the number at www.sensica.com).

Sensica requires proof of purchase. Please provide an invoice or purchase receipt showing the date of purchase of the device



Enter your serial number on our website to activate your warranty.
Saisissez votre numéro de série sur notre site internet pour activer votre garantie.
Inserisci il tuo numero di serie sul nostro sito Web per attivare la garanzia.
Geben Sie zur Garantieregistrierung Ihre Seriennummer auf unserer Website ein.
Introduzca su número de serie en nuestro sitio web para activar su garantía.

www.sensica.com